



NADAL[®] SARS-CoV-2/ Influenza A+B/RSV Ag Test

Antigen rapid test for the detection of SARS-CoV-2,
Influenza A & B and RSV



Sample material	Sensitivity	Specificity	Product code	Pack size
NADAL [®] SARS-CoV-2/Influenza A+B/RSV Ag Test				
Nasopharyngeal swab	Influenza A: 94.83% Influenza B: 96.97% RSV: 98.23% SARS-CoV-2: 96.23%	Influenza A: 98.68% Influenza B: 98.48% RSV: 99.20% SARS-CoV-2: 99.61%	243502N-20	20 tests per pack

Product-Features:

- ✓ Testing of four pathogens simultaneously
- ✓ High sensitivity
- ✓ Fast results in only 15 minutes

Test kit contains

- 20 test cassettes
- 20 sterile swabs
- 20 extraction tubes with buffer and dropper caps
- 1 reagent holder
- 1 package insert



NADAL[®] SARS-CoV-2/ Influenza A+B/RSV Ag Test



4-in-1 Test


A number of viruses can cause infections of the respiratory tract, some of which can be severe. In high-risk patients, it is often important to determine which pathogen is causing the illness. With this new combination rapid test, it is possible to test for SARS-CoV-2, Influenza A & B and RSV from within your practice, meaning no time is wasted waiting for laboratory results.

NADAL[®] SARS-CoV-2/Influenza A+B/RSV Ag


SARS-CoV-2/Influenza A+B/RSV Ag Combo Rapid Test is a rapid chromatographic immunoassay for the qualitative detection of SARS-CoV-2, Influenza A, Influenza B and Respiratory Syncytial Virus (RSV) antigens in nasopharyngeal swab specimens. It is intended for clinical laboratories and healthcare professional use only.

We help you care!

Test procedure



Result interpretation



A detailed description can be found in the package insert.

Dansk distributør/forhandler:

Jepsen HealthCare

Tune Erhvervspark, Tune Parkvej 5, DK 4030 Tune

Telefon 44991432 mail@jepsenhealthcare.dk www.jepsenhealthcare.dk



info@nal-vonminden.com
www.nal-vonminden.com